# DEI® RAINBOW

## HARD BITE REGISTRATION

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Polydimethyl Polymethyl Hydrogen Siloxane  
Polydimethyl Methyl Vinylsiloxane / Silica / Paraffine

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# DEI® RAINBOW

## HARD BITE REGISTRATION

### VINYLPOLYSILOXANE FOR OCCLUSAL REGISTRATIONS

Product for dental use only. Keep away from children. Date of specifications: 10/2010

#### Available packages

Four 50 ml cartridges (COD. RAIN 8)

#### Characteristics

An extremely rigid material (Shore A 90) with fast hardening (only 1 minute in oral cavity). Owing to its high thixotropic behaviour and viscosity it can be distributed directly onto the dental arch. Such kind of occlusal registration can be kept for months without risk of size deformations.

#### Composition

Polydimethyl Polymethyl Hydrogen Siloxane / Polydimethyl Methyl Vinylsiloxane / Silica / Paraffin

#### Way of use

Place the cartridge on the distributor gun. Replace the sealing cap with a cannula for mixing. Spread the product directly onto the lower arch and make the occlusion within 20 seconds from the beginning of the delivery. Let cure in the oral cavity for 60 seconds. After use, let the cannula onto the cartridge. Replace the cannula only at the moment of further use of the product. This practice avoids contamination of base paste with catalysts causing a plug into the cartridge. In this case, the plug must be eliminated before placing a new cannula.

#### Warning

The use of the product requires the distributor gun and the cannulas for mixing (not included in the pack).

#### Technical data

Shore A hardness: ab. 90 - Size deformation after 24 h: 0.3%

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